



Small and Standard Incision Cataract Surgery

CT SPHERIS™ 203 · Technical Data

Former product name: XL Stabi

Monofocal spheric IOL for Standard Incision Surgery in 3-haptic design



Optic Design:	Monofocal, spheric
Material:	Hydrophilic acrylic (28%)
Optic Diameter:	6.0 mm
Total Diameter:	10.5 mm
Haptic Angulation:	10°
Lens Design:	Single-Piece
Incision Size:	2.85 mm
Company Labeled	
A-Constant:	118.0
Diopter Range:	From +8.0 to +32.0 D +8.0 to +26.0 D, 0.5 D increments +26.0 to +32.0 D, 1.0 D increments
ACD:	4.97
Implantation in:	Bag
Injector:	DK7784/Microglide 4220
